



Aktin Chemicals, Inc.

No.52, Yongfeng Road, Hi-tech Zone, Chengdu
610041, P.R.China

Phone:86-28-85159085 Fax:86-28-85152372

website: www.aktinchem.com Email: info@aktinchem.com

CERTIFICATE OF ANALYSIS

Fucosterol

Description:	Fucosterol 98% min by HPLC		
CAS:	17605-67-3	MFW	558.70
MF	C ₂₉ H ₄₈ O	Manufacture Date:	Aug. 23, 2015
Batch:	FU-150823	Quantity:	11.5g

ITEMS	SPECIFICATION	RESULTS	METHODS
Appearance	White Powder	Conform	Visual
Identification	Positive	Positive	TLC
Fucosterol	>98%	98.17%	HPLC
Loss on Drying	≤5%	0.39%	2H 105°C
Ashes	≤0.5%	0.08%	2H 525°C

Storage: store in cool & dry place, keep away from strong light and heat

Stability: 36 months when properly stored.

Packing: 10g/vial

Quality Control Manager: Ashily Zhou Quality Assurance Manager: Will Zhang